


JMS

BYPASS TUBE

CORONARY SHUNT TUBE



JMS BYPASS TUBE

SPECIAL FEATURES

SMOOTH OPERATION

Unique Design
Offset tube provides large area for anastomosis




JMS





Conventional

Attached guidewire for easy insertion




Easy grip tab with hole



Kink resistant for maximum maneuverability




UTURE



HIGH FLOW RATE

Flow rate characteristics table




SAFETY


High biocompatibility Polyamide

Snake in Technique

Insertion Method (Snake-In Technique)



Specifications



① Head Diameter (mm)	Configuration				
	1.25	1.50	1.75	2.00	2.50
② Head Length (mm)	3	3	3	3	3
③ Tie point (mm)	8.5	8.5	9.0	9.0	9.0
④ Total length (mm)	17	17	21	21	22
Inner diameter (mm)	0.5	0.5	0.7	1.0	1.5
Outer diameter (mm)	0.7	0.7	1.0	1.3	1.8

Distributor :

JMS North America Corporation

22320 Foothill Blvd., Suite 350
Hayward, CA 94541 USA

Tel : 510-888-9090
Fax : 510-888-9099

Website : www.jmsna.net

Manufacturer :

JMS Co.,Ltd.

12-17 Kako-Machi, Naka-ku
Hiroshima 730-8652, Japan

Tel : +81-(0) 82-243-5505
Fax : +81-(0) 82-243-1976