

Harmony®

JMS Harmony® AV Fistula needle set is used with button-hole access sites. The integrated cap design provides a method for aseptically preparing the access for cannulation.




www.jmsna.net

JMS

Special Features



Buttonhole Bevel



Cap for aseptically removing the buttonhole scab



Color-coded for easy identification

Ordering Information

Single Set

Item Code	Description	Each / Case
827-1400-33	Harmony® 14GX1" BE 12" W/CLAMP	300 ea / case
827-1500-33	Harmony® 15GX1" BE 12" W/CLAMP	300 ea / case
827-1600-33	Harmony® 16GX1" BE 12" W/CLAMP	300 ea / case
827-1700-33	Harmony® 17GX1" BE 12" W/CLAMP	300 ea / case
827-1800-33	Harmony® 18GX1" BE 12" W/CLAMP	300 ea / case

Item Code	Description	Each / Case
827-1402-33	Harmony® 14GX1" BE 16" W/CLAMP	200 ea / case
827-1503-33	Harmony® 15GX1" BE 16" W/CLAMP	200 ea / case
827-1603-33	Harmony® 16GX1" BE 16" W/CLAMP	200 ea / case

*Other configurations are available upon request

For more information please contact us or your local distributor

Copyright©2010 JMS North America Corporation.
All rights reserved.

BRO/HAR/001_R2

Imported By: JMS North America Corporation

Main Office
22320 Foothill Blvd., Suite 350
Hayward, CA 94541
Tel: 510-888-9090
Fax: 510-888-9099
NA: customerservice@jmsna.net

Sales & Marketing
480 Sawgrass Corporate Pkwy., Suite 120
Sunrise, FL 33325
Tel: 954-689-9280
Fax: 954-689-9289
LatAm: ventas@jmsna.net