


Disposable Vein Stripper

Provides a Minimally Invasive Therapy for
Removal of Varicose Veins



Unique shaped tip easily navigates
valves & tortuous veins

The nylon coated stainless steel wire
allows for safe insertion with minimal
trauma

Mid-wire ligation point assures
successful invagination

Simple and secure connection for
optional Babcock technique



Accessories

Phlebectomy Hook



- Handle provides clear indication of hook position
- Blunt tip prevents accidental trauma to the vessel & skin

Tumescent Local Anesthesia (TLA) Needle



- Side Holes allow for uniform irrigation
- Curved cannula adapts to the patient's anatomy
- Available in 14, 15 & 16G
- Length 10"

Bell-Shape Heads and Wires



- Multiple sizes offer the ability to adapt to the vein size (7mm-9mm & 9mm-12mm)
- Stainless steel wire for optional Babcock Technique

JMS Disposable Vein Stripper Product Line:

CODE	PRODUCT DESCRIPTION	REMARKS
142138025	Vein Stripper 2.0 mm	Includes Bell-Shape Heads and Wires (7mm-9mm)
142138015	Vein Stripper 2.5 mm	Includes Bell-Shape Heads and Wires (9mm-12mm)
511420201	TLA Needle 14G	—
511520201	TLA Needle 15G	—
511620201	TLA Needle 16G	—
145640015	Phlebectomy Hook	Available in one size only

For Single Use Only

Sterilized by ETO

3-year Shelf Life

Blister Packaging



For more information please contact us or your local distributor

Copyright ©2008 JMS North America Corporation. All rights reserved.

BRO/VS-E/001_R1

Imported By:

JMS North America Corporation

MAIN OFFICE

22320 Foothill Blvd., Suite 350
Hayward, CA 94541, U.S.A.

Tel: +1-510-888-9090

Fax: +1-510-888-9099

NA: customerservice@jmsna.net

SALES & MARKETING

480 Sawgrass Corporate Pkwy, Suite 120
Sunrise, FL 33325, U.S.A.

Tel: +1-954-689-9280

Fax: +1-954-689-9289

LatAm: ventas@jmsna.net